

PROPHY PRODUCTS



DTP00730 – VETCARE
Prophy Paste
250 gr. Jar Fruit Flavor

DTP00729 – VETCARE
Prophy Paste
Mint Flavor

DTP00735 – Prophy
Paste Finger Ring Cup
Stainless Steel



New 'TWISTER' DISPOSABLE PROPHY ANGLE

Prophy Cup reciprocates back and forth.
No paste splatter! No heat!
No more tangles from catching hairs!

DTP10967 – Twister Prophy Angle with Cup, Bag of 100
DTP00731 – Individual Prophy Paste Cups, Box of 100
DTP10965 – Standard Disposable Prophy Angle with Cup, box of 144



ORALDENT CLEANSING SOLUTION

DA100CS


0.1% Chlorhexidine solution. 8 oz.



INDIVIDUAL PROPHY PAST CUPS

Package of 100

DTP00731



TRUE-SEAL


Prophy Angle
(Sealed & Autoclavable)

DA10955 – Snap-on

ECONOMY

Prophy Angle
(Autoclavable)

DA551PSN – Snap-on



TRUE-SEAL


Prophy Angle
(Sealed & Autoclavable)

DA10956 – Screw-on

ECONOMY

Prophy Angle
(Autoclavable)

DA552PSC – Screw-on



CLOSE UP

PROPHY CUPS

DTP00720 – Snap-on Prophy Cup, box of 144
DTP00670 – Screw-on Prophy Cup, box of 144
DTP00665 – Latch Type Prophy Cup,




FLUORIDE QUICK-GEL

(Stannous Fluoride .4%)
Helps prevent cavities by
strengthening the enamel.
Also has an anti-bacterial ef-
fect. 4.3 oz. bottle
DTP10522 – Cinnamon
DTP10524 – Berry
DTP10526 – Mint



FLUORIDE CONCENTRATE

1.23% Acidulated &
phosphorilated
DTP10525 – 16 oz. bottle



SAFETY GLASSES

(Clear) DTP360005



SPLASH SHIELD

(Protection from paste
splatter/bacteria)
DTP11234 – Pkg. of 3



DISCLOSING SOLUTION

Used to locate calculus or plaque
on teeth. Good to spot check thor-
oughness of prophylaxis. Paint on
with Q-tip or brush. Plaque Disclosing Solu-
tion
DTP10695 - 2 oz.